



Dosing Innovation for Liquid and Powder Chemicals




#### Main advantage

- Crease-and spot-free sample dyeings
- Knitted and woven fabrics as well as yarn can be dyed
- Natural and synthetic materials
- State of the art process controller with large integrated program memory
- Low liquor ratio: 1:5 fabric dyeing cotton
- Stainless steel cabinet
- Dosing without opening the beakers
- No glycerin or cooling water needed
- Very good price / performance value

## Simple operation

- Filling dye liquor and textile
- Simple fixation of the beakers into the machine
- Calling of the dye program from memory and starting the dye process







# **Beaker support**

- Three-dimensional movement and 360° rotation (clockwise and counterclockwise) allow even and dyed through samples
- The multifunctional beaker holder allows the use of different beaker sizes
- Reversable rotation and speed between 5 and 100 rotations per minute

#### **Beaker sizes**

- 175 ml x 16 beakers (standard: Piece samples only)
- 350 ml x 12 beakers (for yarn or piece samples)
- 500 ml x 8 beakers



# **Dosing Innovation for Liquid and Powder Chemicals**

# Disadvantages of no adding possibility

- Time is lost (opening and closing beaker) and cooling down of the dye bath takes place. This increases the possibility of creases and wrinkles.
- Adding over time (as dosing on the production machines) is not possible.
- Danger of having spots on fabric, because of adding chemicals directly into the dye bath.
- Time consuming operation and great care of the operator is needed. Errors are frequent and reproduction from lab to production is not optimal.

## Disadvantages of the injection method

- Time consumption during adding relatively big, since for every beaker approx. 10 to 15 seconds are being lost (on a 16 beaker machine, adding takes therefore approx. 3 - 4 minutes). This increases the possibility of crease marks and wrinkles on fabrics.
- Only solutions (no powder chemicals) can be added. This increases liquor ratio, e.g. when soda ash has to be added. Reproduction from lab to production is not optimal.
- The solution is given to the bath in a very short time. Spots are very frequent since the chemical is not added in a time frame as on the production machine.
- Membranes and injection device need spare parts and are therefore costly.

# **Our system**

- With the optional specially designed lids, powder as well as liquid chemicals can be put into a small beaker (separated from the dye bath).
- Actuating a lever (time approx. 3 seconds per beaker = less than 1 minute for a 16 beaker machine) initiates the slow dilution of the chemical into the dye bath (beaker lid must not be opened to activate this process).
- Comparing the above systems with our system: Beakers do not move for only one minute. This increases the quality (less wrinkle and crease marks).
- Since the chemical mixes slowly with the dye bath, concentrated chemicals are not touching the fabric and the fabric is dyed spot-free.
- It is the optimal method for the cotton and cotton/polyester dyers but also for any other materials.







#### **Process controller**

- Setex 575 CE or Sedo-Treepoint 1808+ controllers
- Controllers can be connected to Orgatex or Sedomaster
- Graphical display showing dyeing curve, process steps, temperature, speed and rotation direction.

#### **Temperature measuring**

The measuring of the temperature takes place inside the beakers allowing the direct measurement of the dye liquor. This enables the machine to accurately control and regulate the temperature.

## **Heating system**

- The laboratory dye machine is equipped with two 1 kW or 1.5 kW (max.) infrared heating tubes.
- Maximal temperature: 140 °C.
- Heating rate can be varied between 0.1 °C/min and 4°C/min.

# **Cooling system**

Air cools down to the required temperature. The machine does not need any cooling water.





# **Options**

- The standard beaker holder has 16 positions.
  Optionally, a holder with 12 beaker positions is available.
- To dye yarn we suggest our yarn holder. This accessory is available in combination with the 350 ml beakers.
- Teflon foil to dye crease sensitive and delicate materials.

### **Specifications**

810 x 650 x 910 mm
120 kg
2.5 to 3.5 kW
3 kW
400 V, 50 or 60 Hz
40 °C

LAB-PRO GmbH Dorfstrasse 19, CH-4806 Wikon/Switzerland

Tel.:+41 62 745 16 00 info@lab-pro.ch Fax:+41 62 745 16 01 www.lab-pro.ch

