


# Rotator for IRD-III Beakers



## Main advantages:

- Function: to move up to 32 prepared beakers.
- Less uneven dyeings due to constant moving of the dye beakers after the preparation.
- Higher efficiency in the laboratory since samples can be prepared in advance.
- Better order in the laboratory.
- Good price / performance value.
- Stainless steel.



## Easy usage:

- Plugging in of the beakers.
- Push button to start moving the beakers.

## Possible beaker sizes:

- 175 ml beakers
- 350 ml beakers
- 500 ml beakers

## Specification:

Size (W x D x H):	860 x 450 x 630 mm
Machine weight:	approx. 25 - 30 kg
Max. working temperature:	60 °C

LAB 102.0 e