

LD Laboratory Dryer



Main advantages:

- Dryer to dry textile samples such as yarn, loose fibers and fabric (sometimes our laboratory stenter is needed to match lab to production).
- Temperature range between room temperature and 200 °C.
- Fast heat up due to 6 kW heating elements.
- Wide heating chamber.
- Simple operation.
- Solid and compact table design.
- Plug and play (only electrical connection needed).
- Excellent price value ratio.



Specification:

Dimensions (W x D x H):	820 x 630 x 540 mm
Drying space (W x D x H):	400 x 400 x 350 mm
Drying temperature:	up to 200°C
Machine weight:	80 kg
Power consumption:	6 KW
Voltage:	Single phase 230 V, 50 Hz

LAB 106.0 e