

LDP Laboratory Dye Padder



Specification:

LDP1	adopt high quality latex
LDP2	adopt America Hypalon® latex
Padder pressure:	0.1 - 0.8 Mpa (1 - 8 bar)
Padder length:	400 mm
Padder diameter:	200 mm
Padder hardness:	approx. 70° A
Machine weight:	approx. 100 kg
Dimensions (W x D x H):	820 x 600 x 490 mm
Voltage:	380 V, 50 Hz
Power:	750 W
Wattage:	5 A

LAB 111.0 e

LAB-PRO GmbH
Dorfstrasse 19
CH-4806 Wikon / Switzerland

Phone: +41 62 745 16 00
Fax: +41 62 745 16 01

info@lab-pro.ch
www.lab-pro.ch